

F1-CL-B1 試薬棚 中央流し



F1-CL36-B1

Model	Type	Dimensions(mm)			Weight (approx.kg)	Code
		W	D	H		
F1-CL30-B1	L	3000	1500	800 (1790)	435	311-19113-01
F1-CL36-B1		3600			475	311-19114-01
F1-CL42-B1		4200			592	311-19115-01

本体高さ900(1890)mmのタイプも特注にて対応可能です。
900(1890)mm height types are available on custom order.

仕様

天板	サイファス ポリプロピレン製エッジ 30mm厚
本体	角鋼管フレーム 焼付け塗装 両面化粧パーティクルボード
脚部	レベルアジャスター付
流し	2槽両面式ステンレス流し
試薬棚	B1-W ×2または4 本体：両面化粧パーティクルボード
付属品	ワゴン(E1-A・E1-C×各2) 水栓(台付3方口) LED照明
積載荷重	1スパンあたり200kg(等分布荷重)

●流しTW-C15についてはP.239を、試薬棚B1-WについてはP.233をご参照ください。


Spec.

Table top	Cyphus (synthetic resin), P.P. edging, entire thickness 30mm
Base unit	Steel square pipe frame, enameled finish Both sides decorative particle board
Leg part	Level adjusters attached
Sink	One double-sided stainless sink with two drain traps
Reagent shelf	B1-W x2 or 4 Base unit: Both sides decorative particle board
Accessory	Wagon (E1-A, E1-C x2 each) Water faucet (3 way table top type x1) LED lamp
Weight capacity	200kg per 1 span (even distribution weight capacity)


●See P.239 for sink TW-C15, and P.233 for reagent shelf B1-W.

寸法図


F1-CL30-B1




F1-CL36-B1



F1-CL42-B1




Type : L



配管立ち上がり

Service and roughing-in details



※点線内：配管立ち上がり許容範囲
W 給水 D 排水 AC 電気

*In the dotted-line: Piping space
W : Water D : Waste AC : E 100V