

F1-NL-B1 試薬棚



F1-NL30-B1

Model	Type	Dimensions(mm)			Weight (approx.kg)	Code
		W	D	H		
F1-NL18-B1	L	1800	1500	800 (1790)	316	311-18432-01
F1-NL24-B1		2400			381	311-18433-01
F1-NL30-B1		3000			421	311-18434-01
F1-NL36-B1		3600			538	311-18435-01

本体高さ900(1890)mmのタイプも特注にて対応可能です。
900(1890)mm height types are available on custom order.

仕様

天板	サイファス ポリプロピレン製エッジ 30mm厚
本体	角鋼管フレーム 焼付け塗装 両面化粧パーティクルボード
脚部	レベルアジャスター付
試薬棚	B1-W ×2または4 本体：両面化粧パーティクルボード
付属品	ワゴン(E1-A・E1-C×各2) LED照明
積載荷重	1スパンあたり200kg(等分布荷重)


●試薬棚B1-WについてはP.233をご参照ください。

Spec.

Table top	Cyphus (synthetic resin), P.P. edging, entire thickness 30mm
Base unit	Steel square pipe frame, enameled finish Both sides decorative particle board
Leg part	Level adjusters attached
Reagent shelf	B1-W x2 or 4 Base unit: Both sides decorative particle board
Accessory	Wagon (E1-A, E1-C x2 each) LED lamp
Weight capacity	200kg per 1 span (even distribution weight capacity)


●See P.233 for more information about the reagent shelf B1-W.

寸法図



配管立ち上がり

Service and roughing-in details



※点線内：配管立ち上がり許容範囲
AC 電気

*In the dotted-line: Piping space
AC: E 100V