

# F1-SL-SB1 試薬棚 スポット流し ステンレス流し



F1-SL36-SB1

| Model       | Type | Dimensions (mm) |      |            | Weight (approx.kg) | Code         |
|-------------|------|-----------------|------|------------|--------------------|--------------|
|             |      | W               | D    | H          |                    |              |
| F1-SL24-SB1 | L    | 2400            | 1500 | 800 (1790) | 383                | 311-18464-01 |
| F1-SL30-SB1 |      | 3000            |      |            | 452                | 311-18465-01 |
| F1-SL36-SB1 |      | 3600            |      |            | 476                | 311-18466-01 |
| F1-SL42-SB1 |      | 4200            |      |            | 571                | 311-18467-01 |

本体高さ900(1890)mmのタイプも特注にて対応可能です。  
 900(1890)mm height types are available on custom order.

## 仕様

|      |  |
|------|--|
| 天板   | サイファス ポリプロピレン製エッジ 30mm厚                                    |
| 本体   | 角鋼管フレーム 焼付け塗装<br>両面化粧パーティクルボード                             |
| 脚部   | レベルアジャスター付   |
| 流し   | 1槽ステンレス流し<br>スポット流し ×1または2                                 |
| 試薬棚  | B1-W ×2または3<br>本体：両面化粧パーティクルボード                            |
| 付属品  | ワゴン(E1-A・E1-C×各2)<br>水栓(台付3方口)<br>(ゲースネック式×2または4)<br>LED照明 |
| 積載荷重 | 1スパンあたり200kg(等分布荷重)  |


●流しTW-S15についてはP.239を、試薬棚B1-WについてはP.233をご参照ください。

## Spec.

|                 |   |
|-----------------|---|
| Table top       | Cyphus (synthetic resin), P.P. edging, entire thickness 30mm  |
| Base unit       | Steel square pipe frame, enameled finish<br>Both sides decorative particle board                                |
| Leg part        | Level adjusters attached  |
| Sink            | One single-sided stainless sink with a drain trap<br>Stainless steel (SUS-304) spot sink with P.P. trap x1 or 2 |
| Reagent shelf   | B1-W x2 or 3<br>Base unit: Both sides decorative particle board<br>Wagon (E1-A, E1-C x2 each)                   |
| Accessory       | Water faucet (3 way table top type x1)<br>(Gooseneck table top type x2 or 4)<br>LED lamp                        |
| Weight capacity | 200kg per 1 span (even distribution weight capacity)  |


●See P.239 for sink TW-S15, and P.233 for reagent shelf B1-W.

## 寸法図



## 配管立ち上がり

### Service and roughing-in details



※点線内：配管立ち上がり許容範囲  
 W 給水 D 排水 AC 電気  
 \*In the dotted-line: Piping space  
 W: Water D: Waste AC: E 100V

| 開口                 | a   | b    |
|--------------------|-----|------|
| Both sided cabinet |     |      |
| 3000               | 600 | 750  |
| 3600               | 750 | 900  |
| 4200               | 950 | 1150 |