

G1-CL-B1T

試薬棚 試験管掛け 中央流し



G1-CL36-B1T

Model	Type	Dimensions (mm)			Weight (approx. kg)	Code
		W	D	H		
G1-CL30-B1T	L	3000	1500	800 (1790)	436	311-18219-01
G1-CL36-B1T		3600			510	311-18220-01
G1-CL42-B1T		4200			589	311-18221-01

本体高さ900(1890)mmのタイプも特注にて対応可能です。
 900(1890)mm height types are available on custom order.

仕様

天板	サイファス ポリプロピレン製エッジ 30mm厚
本体	両面化粧パーティクルボード
流し	2槽両面式ステンレス流し
試薬棚	B1-W ×2または4 本体：両面化粧パーティクルボード
試験管掛け	両面 TA-W 1台
付属品	水栓(壁付3方口×2) LED照明

●流しTW-C15HについてはP.239を、試薬棚B1-WについてはP.233を、試験管掛けTA-WについてはP.238をご参照ください。


Spec.

Table top	Cyphus (synthetic resin), P.P. edging, entire thickness 30mm
Base unit	Both sides decorative particle board
Sink	One double-sided stainless sink with two drain traps
Reagent shelf	B1-W x2 or 4 Base unit: Both sides decorative particle board
Peg board	TA-W both side type x1
Accessory	Water faucet (3 way wall mounted type x2) LED lamp


●See P.239 for sink TW-C15H, P.233 for reagent shelf B1-W, and P.238 for peg board TA-W.

寸法図


G1-CL30-B1T



G1-CL36-B1T




G1-CL42-B1T




配管立ち上がり

Service and roughing-in details



Type : L



※点線内：配管立ち上がり許容範囲
 W 給水 D 排水 AC 電気
 *In the dotted-line: Piping space
 W: Water D: Waste AC: E 100V