

# G1-CL-B3 試薬棚 中央流し



G1-CL36-B3

Model	Type	Dimensions (mm)			Weight (approx.kg)	Code
		W	D	H		
G1-CL30-B3	L	3000	1500	800 (1790)	409	311-18240-01
G1-CL36-B3		3600			479	311-18241-01
G1-CL42-B3		4200			576	311-18242-01

本体高さ900(1890)mmのタイプも特注にて対応可能です。  
 900(1890)mm height types are available on custom order.

## 仕様

天板	サイファス ポリプロピレン製エッジ 30mm厚
本体	両面化粧パーティクルボード
流し	2槽両面式ステンレス流し フラットタイプ
試薬棚	B3-W ×2または4 本体：両面化粧パーティクルボード
付属品	水栓(台付2方口×2)
積載荷重	100kg/m <sup>2</sup>

●試薬棚B3-WについてはP.234をご参照ください。


## Spec.

Table top	Cyphus (synthetic resin), P.P. edging, entire thickness 30mm
Base unit	Both sides decorative particle board
Sink	One double-sided stainless sink with two drain traps
Reagent shelf	B3-W x2 or 4 Base unit: Both sides decorative particle board
Accessory	Water faucet (2 way table top type x2)
Weight capacity	100kg/m <sup>2</sup>

●See P.234 for more information about the reagent shelf B3-W.

## 寸法図


G1-CL30-B3



G1-CL36-B3




G1-CL42-B3



## 配管立ち上がり

Service and roughing-in details



※点線内：配管立ち上がり許容範囲  
 W 給水 D 排水 AC 電気

\*In the dotted-line: Piping space  
 W: Water D: Waste AC: E 100V

## Type: L

