

# G1-CL-B5 試薬棚 中央流し



G1-CL36-B5

Model	Type	Dimensions(mm)			Weight (approx.kg)	Code
		W	D	H		
G1-CL30-B5	L	3000	1500	800 (1790)	425	311-18270-01
G1-CL36-B5		3600			497	311-18271-01
G1-CL42-B5		4200			604	311-18272-01

本体高さ900(1890)mmのタイプも特注にて対応可能です。  
900(1890)mm height types are available on custom order.

## 仕様

天板	サイファス ポリプロピレン製エッジ 30mm厚
本体	両面化粧パーティクルボード
流し	2槽両面式ステンレス流し フラットタイプ
試薬棚	B5-W ×2または4 本体：両面化粧パーティクルボード
付属品	水栓(台付2方口×2) LED照明
積載荷重	100kg/m <sup>2</sup>


●試薬棚B5-WについてはP.235をご参照ください。

## Spec.

Table top	Cyphus (synthetic resin), P.P. edging, entire thickness 30mm
Base unit	Both sides decorative particle board
Sink	One double-sided stainless sink with two drain traps
Reagent shelf	B5-W x2 or 4 Base unit: Both sides decorative particle board
Accessory	Water faucet (2 way table top type x2) LED lamp
Weight capacity	100kg/m <sup>2</sup>


●See P.235 for more information about the reagent shelf B5-W.

## 寸法図



## 配管立ち上がり

Service and roughing-in details



※点線内：配管立ち上がり許容範囲  
W 給水 D 排水 AC 電気

\*In the dotted-line: Piping space  
W: Water D: Waste AC: E 100V

