

G1-NL-SB5 試薬棚 スポット流し



G1-NL30-SB5

Model	Type	Dimensions (mm)			Weight (approx.kg)	Code
		W	D	H		
G1-NL18-SB5	L	1800	1500	800	305	311-19106-01
G1-NL24-SB5		2400		(1790)	367	311-19107-01
G1-NL30-SB5		3000		447	311-19108-01	

本体高さ900(1890)mmのタイプも特注にて対応可能です。
 900(1890)mm height types are available on custom order.

仕様

天板	サイファス ポリプロピレン製エッジ 30mm厚
本体	両面化粧パーティクルボード
流し	スポット流し ×2
試薬棚	B5-W ×2 本体：両面化粧パーティクルボード
付属品	水栓(直水栓×4) LED照明

●試薬棚B5-WについてはP.235をご参照ください。

Spec.

Table top	Cyphus (synthetic resin), P.P. edging, entire thickness 30mm
Base unit	Both sides decorative particle board
Sink	Stainless steel (SUS-304) spot sink with P.P. trap x2
Reagent shelf	B5-W x2 Base unit: Both sides decorative particle board
Accessory	Water faucet (Direct water faucet x4) LED lamp


●See P.235 for more information about the reagent shelf B5-W.

寸法図


G1-NL18-SB5




G1-NL24-SB5



G1-NL30-SB5




Type : L



配管立ち上がり

Service and roughing-in details



※点線内：配管立ち上がり許容範囲
 W 給水 D 排水 AC 電気

*In the dotted-line : Piping space
 W : Water D : Waste AC : E 100V

間口	a
Both sided cabinet	
1800	300
2400	600
3000	750