

G1-SL-B1

試薬棚

ステンレス流し



G1-SL36-B1

概要・ガイド

スチール製実験台

中央実験台

フレーム中央実験台

サイド実験台

フレームサイド実験台

実験台用ユーティリティ

装置台・昇降台・作業台

除振台・天秤台

流し台

関連オプション

01

ドラフトチャンパー

02

排ガス処理装置

03

排風機

04

実験台

05

戸棚・薬品庫・ワゴン

06

クリーン・バスタード・環境試験施設/機器

07

動物関連施設/機器

08

医学研究・検査施設

09

廃水処理装置

10

機器移設・メンテナンス

11

技術資料

12

さくいん

Model	Type	Dimensions (mm)			Weight (approx.kg)	Code
		W	D	H		
G1-SL24-B1	L	2400	1500	800 (1790)	354	311-18194-01
G1-SL30-B1		3000			416	311-18195-01
G1-SL36-B1		3600			496	311-18196-01
G1-SL42-B1		4200			575	311-18197-01

本体高さ900(1890)mmのタイプも特注にて対応可能です。
900(1890)mm height types are available on custom order.

仕様

天板	サイファス ポリプロピレン製エッジ 30mm厚
本体	両面化粧パーティクルボード
流し	1槽ステンレス流し
試薬棚	B1-W ×2または3 本体：両面化粧パーティクルボード
付属品	水栓(台付3方口) LED照明
積載荷重	100kg/m ²

●試薬棚B1-WについてはP.233をご参照ください。


Spec.

Table top	Cyphus (synthetic resin), P.P. edging, entire thickness 30mm
Base unit	Both sides decorative particle board
Sink	One single-sided stainless sink with a drain trap
Reagent shelf	B1-W x2 or 3 Base unit: Both sides decorative particle board
Accessory	Water faucet (3 way table top type x1) LED lamp
Weight capacity	100kg/m ²


●See P.233 for more information about the reagent shelf B1-W.

寸法図


G1-SL24-B1




G1-SL36-B1




G1-SL30-B1



G1-SL42-B1




Type : L



配管立ち上がり

Service and roughing-in details



※点線内：配管立ち上がり許容範囲
W 給水 D 排水 AC 電気
*In the dotted-line: Piping space
W: Water D: Waste AC: E 100V