

G1-TL-B5 試薬棚 陶器流し



G1-TL30-B5

Model	Type	Dimensions(mm)			Weight (approx.kg)	Code
		W	D	H		
G1-TL18-B5	L	1800	1500	800 (1790)	322	311-18262-01
G1-TL24-B5		2400			384	311-18263-01
G1-TL30-B5		3000			464	311-18264-01
G1-TL36-B5		3600			546	311-18265-01

本体高さ900(1890)mmのタイプも特注にて対応可能です。
900(1890)mm height types are available on custom order.

仕様

天板	サイファス ポリプロピレン製エッジ	30mm厚
本体	両面化粧パーティクルボード	
流し	陶器流し	
試薬棚	B5-W ×2または3 本体：両面化粧パーティクルボード	
付属品	水栓(台付3方口) LED照明	
積載荷重	100kg/m ²	

●試薬棚B5-WについてはP.235をご参照ください。


Spec.

Table top	Cyphus (synthetic resin), P.P. edging, entire thickness 30mm
Base unit	Both sides decorative particle board
Sink	Ceramic sink B5-W x2 or 3
Reagent shelf	Base unit: Both sides decorative particle board
Accessory	Water faucet (3 way table top type x1) LED lamp
Weight capacity	100kg/m ²


●See P.235 for more information about the reagent shelf B5-W.

寸法図


G1-TL18-B5




G1-TL30-B5




G1-TL24-B5



G1-TL36-B5




Type : L



配管立ち上がり

Service and roughing-in details



※点線内：配管立ち上がり許容範囲
W 給水 D 排水 AC 電気

*In the dotted-line: Piping space
W : Water D : Waste AC : E 100V