

G1-TL-SB5

試薬棚

スポット流し

陶器流し



G1-TL30-SB5

Model	Type	Dimensions(mm)			Weight (approx.kg)	Code
		W	D	H		
G1-TL18-SB5	L	1800	1500	800 (1790)	330	311-18302-01
G1-TL24-SB5		2400			392	311-18303-01
G1-TL30-SB5		3000			472	311-18304-01

本体高さ900(1890)mmのタイプも特注にて対応可能です。
900(1890)mm height types are available on custom order.

仕様

天板	サイファス ポリプロピレン製エッジ 30mm厚
本体	両面化粧パーティクルボード
流し	陶器流し スポット流し ×2
試薬棚	B5-W ×2 本体：両面化粧パーティクルボード
付属品	水栓(台付3方口) (直水栓×4) LED照明

●試薬棚B5-WについてはP.235をご参照ください。


Spec.

Table top	Cyphus (synthetic resin), P.P. edging, entire thickness 30mm
Base unit	Both sides decorative particle board
Sink	Ceramic sink x1 Stainless steel (SUS-304) spot sink with P.P. trap x2
Reagent shelf	B5-W x2 Base unit: Both sides decorative particle board
Accessory	Water faucet (3 way table top type x1) (Direct water faucet x4) LED lamp


●See P.235 for more information about the reagent shelf B5-W.

寸法図


G1-TL18-SB5




G1-TL24-SB5



G1-TL30-SB5




Type : L



配管立ち上がり

Service and roughing-in details



※点線内：配管立ち上がり許容範囲
W 給水 D 排水 AC 電気

*In the dotted-line : Piping space
W : Water D : Waste AC : E 100V

開口	
Both sided cabinet	a
1800	300
2400	600
3000	750